

Antibodies Product Datasheet

Product Name:	LWS Antibody
Catalog Number:	C41081
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LWS. at AA range: 380-420
Concentration:	1 mg/ml
Gene Name:	RXRG NR2B3
Protein Name:	LWS
Other Name:	Retinoic acid receptor RXR-gamma (Nuclear receptor subfamily 2 group B member 3) (Retinoid X receptor gamma)
Human Swiss Prot No.:	P48443
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LWS
Storage and Stability:	-20°C/1 year