

Antibodies Product Datasheet

Product Name:	MT-ATP8 Antibody
Catalog Number:	C41137
Reactivity:	Human, Rat, Mouse
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MT-ATP8
Concentration:	1 mg/ml
Gene Name:	MT-ATP8 ATP8 ATPASE8 MTATP8
Protein Name:	MT-ATP8
Other Name:	ATP synthase protein 8 (A6L) (F-ATPase subunit 8)
Human Swiss Prot No.:	P03928
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human MT-ATP8
Storage and Stability:	-20°C/1 year