

Antibodies Product Datasheet

Product Name:	NDUFA1 Antibody
Catalog Number:	C41169
Reactivity:	Human, Mouse, Rat
Applications:	IHC, ELISA
Dilutions:	IHC-p 1:50-200, ELISA(peptide)1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NDUFA1
Concentration:	1 mg/ml
Gene Name:	NDUFA1
Protein Name:	NDUFA1
Other Name:	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 1 (Complex I-MWFE) (CI-MWFE) (NADH-ubiquinone oxidoreductase MWFE subunit)
Human Swiss Prot No.:	O15239
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human NDUFA1
Storage and Stability:	-20°C/1 year