

Antibodies Product Datasheet

Product Name:	Rabbit anti alpaca IgG Antibody
Catalog Number:	C41396
Reactivity:	Alpaca
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from Alpaca IgG
Concentration:	1 mg/ml
Gene Name:	Rabbit anti alpaca IgG
Protein Name:	Rabbit anti alpaca IgG
Other Name:	
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Alpaca Rabbit anti alpaca IgG
Storage and Stability:	-20°C/1 year