

Antibodies Product Datasheet

Product Name:	CAC1E Antibody
Catalog Number:	C41732
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 370-450
Concentration:	1 mg/ml
Gene Name:	CACNA1E CACH6 CACNL1A6
Protein Name:	Voltage-dependent R-type calcium channel subunit alpha-1E (Brain calcium channel II) (BII) (Calcium channel, L type, alpha-1 polypeptide, isoform 6) (Voltage-gated calcium channel subunit alpha Cav2.3)
Other Name:	
Human Swiss Prot No.:	Q15878
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CAC1E Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year