

Antibodies Product Datasheet

Product Name:	CHD8 Antibody
Catalog Number:	C41735
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CHD8 HELSNF1 KIAA1564
Protein Name:	Chromodomain-helicase-DNA-binding protein 8 (CHD-8) (EC 3.6.4.12) (ATP-dependent helicase CHD8) (Helicase with SNF2 domain 1)
Other Name:	
Human Swiss Prot No.:	Q9HCK8
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CHD8 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year