

Antibodies Product Datasheet

Product Name:	CP250 Antibody
Catalog Number:	C41740
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CEP250 CEP2 CNAP1
Protein Name:	Centrosome-associated protein CEP250 (250 kDa centrosomal protein) (Cep250) (Centrosomal Nek2-associated protein 1) (C-Nap1) (Centrosomal protein 2)
Other Name:	
Human Swiss Prot No.:	Q9BV73
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CP250 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year