

Antibodies Product Datasheet

Product Name:	CT47A Antibody
Catalog Number:	C41741
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 210-290
Concentration:	1 mg/ml
Gene Name:	CT47A1 CT47.1; CT47A2 CT47.2; CT47A3 CT47.3; CT47A4 CT47.4; CT47A5 CT47.5; CT47A6 CT47.6; CT47A7 CT47.7; CT47A8 CT47.8; CT47A9 CT47.9; CT47A10 CT47.10; CT47A11 CT47.11; CT47A12 CT47.12
Protein Name:	Cancer/testis antigen 47A (Cancer/testis antigen 47) (CT47)
Other Name:	
Human Swiss Prot No.:	Q5JQC4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CT47A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year