

Antibodies Product Datasheet

Product Name:	GA45B Antibody
Catalog Number:	C41749
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 90-170
Concentration:	1 mg/ml
Gene Name:	GADD45B MYD118
Protein Name:	Growth arrest and DNA damage-inducible protein GADD45 beta (Myeloid differentiation primary response protein MyD118) (Negative growth regulatory protein MyD118)
Other Name:	
Human Swiss Prot No.:	O75293
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GA45B Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year