

Antibodies Product Datasheet

Product Name:	GCP2 Antibody
Catalog Number:	C41750
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	TUBGCP2 GCP2
Protein Name:	Gamma-tubulin complex component 2 (GCP-2) (hGCP2) (Gamma-ring complex protein 103 kDa) (h103p) (hGrip103) (Spindle pole body protein Spc97 homolog) (hSpc97)
Other Name:	
Human Swiss Prot No.:	Q9BSJ2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GCP2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year