

Antibodies Product Datasheet

Product Name:	GORS1 Antibody
Catalog Number:	C41751
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	GORASP1 GOLPH5 GRASP65
Protein Name:	Golgi reassembly-stacking protein 1 (Golgi peripheral membrane protein p65) (Golgi phosphoprotein 5) (GOLPH5) (Golgi reassembly-stacking protein of 65 kDa) (GRASP65)
Other Name:	
Human Swiss Prot No.:	Q9BQQ3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	GORS1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year