

Antibodies Product Datasheet

Product Name:	IF6 Antibody
Catalog Number:	C41756
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 180-260
Concentration:	1 mg/ml
Gene Name:	EIF6 EIF3A ITGB4BP OK/SW-cl.27
Protein Name:	Eukaryotic translation initiation factor 6 (eIF-6) (B(2)GCN homolog) (B4 integrin interactor) (CAB) (p27(BBP))
Other Name:	
Human Swiss Prot No.:	P56537
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	IF6 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year