

Antibodies Product Datasheet

Product Name:	MUC5A Antibody
Catalog Number:	C41763
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 4400-4480
Concentration:	1 mg/ml
Gene Name:	MUC5AC MUC5
Protein Name:	Mucin-5AC (MUC-5AC) (Gastric mucin) (Lewis B blood group antigen) (LeB) (Major airway glycoprotein) (Mucin-5 subtype AC, tracheobronchial) (Tracheobronchial mucin) (TBM) (Fragments)
Other Name:	
Human Swiss Prot No.:	P98088
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MUC5A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year