

Antibodies Product Datasheet

Product Name:	CENPU Antibody
Catalog Number:	C41997
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 160-240
Concentration:	1 mg/ml
Gene Name:	MLF1IP CENPU ICEN24 KLIP1 PBIP1
Protein Name:	Centromere protein U (CENP-U) (Interphase centromere complex protein 24) (KSHV latent nuclear antigen-interacting protein 1) (MLF1-interacting protein) (Polo-box-interacting protein 1) (centromere protein of 50 kDa) (CENP-50)
Other Name:	
Human Swiss Prot No.:	Q71F23
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CENPU Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year