

Antibodies Product Datasheet

Product Name:	ANM5 Antibody
Catalog Number:	C42005
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	PRMT5 HRMT1L5 IBP72 JBP1 SKB1
Protein Name:	Protein arginine N-methyltransferase 5 (EC 2.1.1.-) (72 kDa ICln-binding protein) (Histone-arginine N-methyltransferase PRMT5) (EC 2.1.1.125) (Jak-binding protein 1) (Shk1 kinase-binding protein 1 homolog) (SKB1 homolog) (SKB1Hs) [Cleaved into: Protein ar
Other Name:	
Human Swiss Prot No.:	O14744
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ANM5 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year