

Antibodies Product Datasheet

Product Name:	CCR9 Antibody
Catalog Number:	C42009
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 140-220
Concentration:	1 mg/ml
Gene Name:	CCR9 GPR28
Protein Name:	C-C chemokine receptor type 9 (C-C CKR-9) (CC-CKR-9) (CCR-9) (G-protein coupled receptor 28) (GPR-9-6) (CD antigen CDw199)
Other Name:	
Human Swiss Prot No.:	P51686
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	CCR9 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year