

Antibodies Product Datasheet

Product Name:	CPT1A Antibody
Catalog Number:	C42010
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CPT1A
Concentration:	1 mg/ml
Gene Name:	CPT1A CPT1
Protein Name:	CPT1A
Other Name:	Carnitine O-palmitoyltransferase 1, liver isoform (CPT1-L) (EC 2.3.1.21) (Carnitine O-palmitoyltransferase I, liver isoform) (CPT I) (CPTI-L) (Carnitine palmitoyltransferase 1A)
Human Swiss Prot No.:	P50416
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.346% sodium azide.
Specificity:	This antibody detects endogenous levels of CPT1A at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year