

Antibodies Product Datasheet

Product Name:	DNM3A Antibody
Catalog Number:	C42012
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	DNMT3A
Protein Name:	DNA (cytosine-5)-methyltransferase 3A (Dnmt3a) (EC 2.1.1.37) (DNA methyltransferase HsaIIIA) (DNA MTase HsaIIIA) (M.HsaIIIA)
Other Name:	
Human Swiss Prot No.:	Q9Y6K1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	DNM3A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year