

Antibodies Product Datasheet

Product Name:	HAVR2 Antibody
Catalog Number:	C42015
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	HAVCR2 TIM3 TIMD3
Protein Name:	Hepatitis A virus cellular receptor 2 (HAVcr-2) (T-cell immunoglobulin and mucin domain-containing protein 3) (TIMD-3) (T-cell membrane protein 3) (TIM-3)
Other Name:	
Human Swiss Prot No.:	Q8TDQ0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HAVR2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year