

Antibodies Product Datasheet

Product Name:	ACAD8 Antibody
Catalog Number:	C42050
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ACAD8
Concentration:	1 mg/ml
Gene Name:	ACAD8 ARC42 IBD
Protein Name:	Isobutyryl-CoA dehydrogenase, mitochondrial (EC 1.3.99.-) (Activator-recruited cofactor 42 kDa component) (ARC42) (Acyl-CoA dehydrogenase family member 8) (ACAD-8)
Other Name:	
Human Swiss Prot No.:	Q9UKU7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ACAD8
Storage and Stability:	-20°C/1 year