

Antibodies Product Datasheet

Product Name:	Actin beta Antibody(Caenorhabditis elegans)
Catalog Number:	C42052
Reactivity:	C.elegans
Applications:	WB
Dilutions:	WB 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Purified recombinant Full length Actin beta (Caenorhabditis elegans)
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	Actin beta
Other Name:	Beta Actin,βActin; YM33033
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Actin beta (Caenorhabditis elegans)
Storage and Stability:	-20°C/1 year