

Antibodies Product Datasheet

Product Name:	Actin beta Antibody(Yeast)
Catalog Number:	C42054
Reactivity:	Yeast
Applications:	WB
Dilutions:	WB 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Purified recombinant Full length Actin beta (Caenorhabditis elegans)
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	Actin beta
Other Name:	Beta Actin, β Actin; YM33042
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Actin beta (Caenorhabditis elegans)
Storage and Stability:	-20°C/1 year