

Antibodies Product Datasheet

Product Name:	ADRO Antibody
Catalog Number:	C42057
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ADRO
Concentration:	1 mg/ml
Gene Name:	FDXR ADXR
Protein Name:	NADPH:adrenodoxin oxidoreductase, mitochondrial (AR) (Adrenodoxin reductase) (EC 1.18.1.2) (Ferredoxin--NADP(+) reductase) (Ferredoxin reductase)
Other Name:	
Human Swiss Prot No.:	P22570
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ADRO
Storage and Stability:	-20°C/1 year