

Antibodies Product Datasheet

Product Name:	AGGF1 Antibody
Catalog Number:	C42058
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human AGGF1
Concentration:	1 mg/ml
Gene Name:	AGGF1 GPATC7 GPATCH7 VG5Q
Protein Name:	Angiogenic factor with G patch and FHA domains 1 (Angiogenic factor VG5Q) (hVG5Q) (G patch domain-containing protein 7) (Vasculogenesis gene on 5q protein)
Other Name:	
Human Swiss Prot No.:	Q8N302
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human AGGF1
Storage and Stability:	-20°C/1 year