

Antibodies Product Datasheet

Product Name:	AP3B1 Antibody
Catalog Number:	C42069
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human AP3B1
Concentration:	1 mg/ml
Gene Name:	AP3B1 ADTB3A
Protein Name:	AP-3 complex subunit beta-1 (Adapter-related protein complex 3 subunit beta-1) (Adaptor protein complex AP-3 subunit beta-1) (Beta-3A-adaptin) (Clathrin assembly protein complex 3 beta-1 large chain)
Other Name:	
Human Swiss Prot No.:	O00203
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human AP3B1
Storage and Stability:	-20°C/1 year