

Antibodies Product Datasheet

Product Name:	AP4M1 Antibody
Catalog Number:	C42070
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human AP4M1
Concentration:	1 mg/ml
Gene Name:	AP4M1 MUARP2
Protein Name:	AP-4 complex subunit mu-1 (AP-4 adapter complex mu subunit) (Adapter-related protein complex 4 mu-1 subunit) (Mu subunit of AP-4) (Mu-adaptin-related protein 2) (mu-ARP2) (Mu4-adaptin) (mu4)
Other Name:	
Human Swiss Prot No.:	O00189
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human AP4M1
Storage and Stability:	-20°C/1 year