

Antibodies Product Datasheet

Product Name:	ARK72 Antibody
Catalog Number:	C42073
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ARK72
Concentration:	1 mg/ml
Gene Name:	AKR7A2 AFAR AFAR1 AKR7
Protein Name:	Aflatoxin B1 aldehyde reductase member 2 (EC 1.1.1.n11) (AFB1 aldehyde reductase 1) (AFB1-AR 1) (Aldoketoreductase 7) (Succinic semialdehyde reductase) (SSA reductase)
Other Name:	
Human Swiss Prot No.:	O43488
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ARK72
Storage and Stability:	-20°C/1 year