

Antibodies Product Datasheet

Product Name:	ASNA Antibody
Catalog Number:	C42078
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ASNA
Concentration:	1 mg/ml
Gene Name:	ASNA1 ARSA TRC40
Protein Name:	ATPase ASNA1 (EC 3.6.-.-) (Arsenical pump-driving ATPase) (Arsenite-stimulated ATPase) (Transmembrane domain recognition complex 40 kDa ATPase subunit) (hARSA-I) (hASNA-I)
Other Name:	
Human Swiss Prot No.:	O43681
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ASNA
Storage and Stability:	-20°C/1 year