

Antibodies Product Datasheet

Product Name:	ATG13 Rabbit Antibody
Catalog Number:	C42081
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of ATG13
Concentration:	1 mg/ml
Gene Name:	ATG13 KIAA0652
Protein Name:	Autophagy-related protein 13
Other Name:	Autophagy-related protein 13; YM3783
Human Swiss Prot No.:	O75143
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous ATG13 protein
Storage and Stability:	-20°C/1 year