

Antibodies Product Datasheet

Product Name:	ATG14L Rabbit Antibody
Catalog Number:	C42083
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of ATG14L
Concentration:	1 mg/ml
Gene Name:	ATG14 KIAA0831
Protein Name:	Beclin 1-associated autophagy-related key regulator (Barkor) (Autophagy-related protein 14-like protein) (Atg14L)
Other Name:	Beclin 1-associated autophagy-related key regulator (Barkor) (Autophagy-related protein 14-like protein) (Atg14L); YM3784
Human Swiss Prot No.:	Q6ZNE5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous ATG14L protein
Storage and Stability:	-20°C/1 year