

Antibodies Product Datasheet

Product Name:	ATLA1 Antibody
Catalog Number:	C42087
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ATLA1
Concentration:	1 mg/ml
Gene Name:	ATL1 GBP3 SPG3A
Protein Name:	Atlastin-1 (EC 3.6.5.-) (Brain-specific GTP-binding protein) (GTP-binding protein 3) (GBP-3) (hGBP3) (Guanine nucleotide-binding protein 3) (Spastic paraplegia 3 protein A)
Other Name:	
Human Swiss Prot No.:	Q8WXF7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ATLA1
Storage and Stability:	-20°C/1 year