

Antibodies Product Datasheet

Product Name:	BABA1 Antibody
Catalog Number:	C42090
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BABA1
Concentration:	1 mg/ml
Gene Name:	BABAM1 C19orf62 MERIT40 NBA1 HSPC142
Protein Name:	BRISC and BRCA1-A complex member 1 (Mediator of RAP80 interactions and targeting subunit of 40 kDa) (New component of the BRCA1-A complex)
Other Name:	
Human Swiss Prot No.:	Q9NWW8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human BABA1
Storage and Stability:	-20°C/1 year