

Antibodies Product Datasheet

Product Name:	BEST3 Antibody
Catalog Number:	C42093
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BEST3
Concentration:	1 mg/ml
Gene Name:	BEST3 VMD2L3
Protein Name:	Bestrophin-3 (Vitelliform macular dystrophy 2-like protein 3)
Other Name:	
Human Swiss Prot No.:	Q8N1M1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human BEST3
Storage and Stability:	-20°C/1 year