

Antibodies Product Datasheet

Product Name:	Bmi-1 Rabbit Antibody(F270)
Catalog Number:	C42097
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein derived from human Bmi-1
Concentration:	1 mg/ml
Gene Name:	BMI1 PCGF4 RNF51
Protein Name:	Bmi-1
Other Name:	Polycomb complex protein BMI-1 (Polycomb group RING finger protein 4) (RING finger protein 51); YM3799
Human Swiss Prot No.:	P35226
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Bmi-1
Storage and Stability:	-20°C/1 year