

Antibodies Product Datasheet

Product Name:	Bovine Serum Albumin Antibody
Catalog Number:	C42100
Reactivity:	Cow
Applications:	WB, ELISA
Dilutions:	WB: 1:2000-5000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Bovine Serum Albumin
Concentration:	
Gene Name:	
Protein Name:	
Other Name:	YM3253
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous BSA protein.
Storage and Stability:	-20°C/1 year