

Antibodies Product Datasheet

Product Name:	C1AS1 Antibody
Catalog Number:	C42104
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human C1AS1
Concentration:	1 mg/ml
Gene Name:	C17orf76-AS1 C17orf45 NCRNA00188
Protein Name:	Putative uncharacterized protein C17orf76-AS1, mitochondrial (C17orf76 antisense RNA 1) (C17orf76 antisense gene protein 1)
Other Name:	
Human Swiss Prot No.:	Q8N1F1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human C1AS1
Storage and Stability:	-20°C/1 year