

Antibodies Product Datasheet

Product Name:	CAR17 Antibody
Catalog Number:	C42107
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CAR17
Concentration:	1 mg/ml
Gene Name:	CARD17 INCA
Protein Name:	Caspase recruitment domain-containing protein 17 (Caspase-1 inhibitor INCA) (Inhibitory caspase recruitment domain protein)
Other Name:	
Human Swiss Prot No.:	Q5XLA6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CAR17
Storage and Stability:	-20°C/1 year