

Antibodies Product Datasheet

Product Name:	CLCA1 Antibody
Catalog Number:	C42126
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CLCA1
Concentration:	1 mg/ml
Gene Name:	CLCA1 CACC1
Protein Name:	Calcium-activated chloride channel regulator 1 (Calcium-activated chloride channel family member 1) (hCLCA1) (Calcium-activated chloride channel protein 1) (CaCC-1) (hCaCC-1)
Other Name:	
Human Swiss Prot No.:	A8K7I4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CLCA1
Storage and Stability:	-20°C/1 year