

Antibodies Product Datasheet

Product Name:	CQ085 Antibody
Catalog Number:	C42140
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CQ085
Concentration:	1 mg/ml
Gene Name:	C17orf85
Protein Name:	Uncharacterized protein C17orf85 (Protein ELG)
Other Name:	
Human Swiss Prot No.:	Q53F19
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CQ085
Storage and Stability:	-20°C/1 year