

Antibodies Product Datasheet

Product Name:	CT47B Antibody
Catalog Number:	C42148
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CT47B
Concentration:	1 mg/ml
Gene Name:	CT47B1 CT47.13
Protein Name:	Cancer/testis antigen 47B (Cancer/testis antigen 47.13) (CT47.13)
Other Name:	
Human Swiss Prot No.:	P0C2W7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CT47B
Storage and Stability:	-20°C/1 year