

Antibodies Product Datasheet

Product Name:	CTRB1 Antibody
Catalog Number:	C42150
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CTRB1
Concentration:	1 mg/ml
Gene Name:	CTRB1 CTRB
Protein Name:	Chymotrypsinogen B (EC 3.4.21.1) [Cleaved into: Chymotrypsin B chain A; Chymotrypsin B chain B; Chymotrypsin B chain C]
Other Name:	
Human Swiss Prot No.:	P17538
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human CTRB1
Storage and Stability:	-20°C/1 year