

Antibodies Product Datasheet

Product Name:	D2HDH Antibody
Catalog Number:	C42155
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human D2HDH
Concentration:	1 mg/ml
Gene Name:	D2HGDH D2HGD
Protein Name:	D-2-hydroxyglutarate dehydrogenase, mitochondrial (EC 1.1.99.-)
Other Name:	
Human Swiss Prot No.:	Q8N465
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human D2HDH
Storage and Stability:	-20°C/1 year