

Antibodies Product Datasheet

Product Name:	DB127 Antibody
Catalog Number:	C42157
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DB127
Concentration:	1 mg/ml
Gene Name:	DEFB127 C20orf73 DEFB27 UNQ1956/PRO6071
Protein Name:	Beta-defensin 127 (Beta-defensin 27) (DEFB-27) (Defensin, beta 127)
Other Name:	
Human Swiss Prot No.:	Q9H1M4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DB127
Storage and Stability:	-20°C/1 year