

Antibodies Product Datasheet

Product Name:	DUS11 Antibody
Catalog Number:	C42170
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DUS11
Concentration:	1 mg/ml
Gene Name:	DUSP11 PIR1
Protein Name:	RNA/RNP complex-1-interacting phosphatase (EC 3.1.3.-) (Dual specificity protein phosphatase 11) (Phosphatase that interacts with RNA/RNP complex 1)
Other Name:	
Human Swiss Prot No.:	O75319
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DUS11
Storage and Stability:	-20°C/1 year