

Antibodies Product Datasheet

Product Name:	DUS8 Antibody
Catalog Number:	C42171
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human DUS8
Concentration:	1 mg/ml
Gene Name:	DUSP8 C11orf81 VH5
Protein Name:	Dual specificity protein phosphatase 8 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase hVH-5)
Other Name:	
Human Swiss Prot No.:	Q13202
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human DUS8
Storage and Stability:	-20°C/1 year