

Antibodies Product Datasheet

Product Name:	EI2BA Antibody
Catalog Number:	C42177
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EI2BA
Concentration:	1 mg/ml
Gene Name:	EIF2B1 EIF2BA
Protein Name:	Translation initiation factor eIF-2B subunit alpha (eIF-2B GDP-GTP exchange factor subunit alpha)
Other Name:	
Human Swiss Prot No.:	Q14232
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human EI2BA
Storage and Stability:	-20°C/1 year