

Antibodies Product Datasheet

Product Name:	EIF3M Antibody
Catalog Number:	C42179
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EIF3M
Concentration:	1 mg/ml
Gene Name:	EIF3M HFLB5 PCID1 GA17 PNAS-125
Protein Name:	Eukaryotic translation initiation factor 3 subunit M (eIF3m) (Fetal lung protein B5) (hFL-B5) (PCI domain-containing protein 1)
Other Name:	
Human Swiss Prot No.:	Q7L2H7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human EIF3M
Storage and Stability:	-20°C/1 year