

Antibodies Product Datasheet

Product Name:	EMAL6 Antibody
Catalog Number:	C42182
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EMAL6
Concentration:	1 mg/ml
Gene Name:	EML6 EML5L
Protein Name:	Echinoderm microtubule-associated protein-like 6 (EMAP-6) (Echinoderm microtubule-associated protein-like 5-like)
Other Name:	
Human Swiss Prot No.:	Q6ZMW3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human EMAL6
Storage and Stability:	-20°C/1 year