

Antibodies Product Datasheet

Product Name:	ENW1 Antibody
Catalog Number:	C42186
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ENW1
Concentration:	1 mg/ml
Gene Name:	ERVW-1 ERVWE1
Protein Name:	HERV-W_7q21.2 provirus ancestral Env polyprotein (Endogenous retrovirus group W member 1) (Env-W) (Envelope polyprotein gp73) (Enverin) (HERV-7q Envelope protein) (HERV-W envelope protein) (Syncytin) (Syncytin-1) [Cleaved into: Surface protein (SU) (gp50
Other Name:	
Human Swiss Prot No.:	Q9UQF0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ENW1
Storage and Stability:	-20°C/1 year