

Antibodies Product Datasheet

Product Name:	ERO1B Antibody
Catalog Number:	C42187
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ERO1B
Concentration:	1 mg/ml
Gene Name:	ERO1LB
Protein Name:	ERO1-like protein beta (ERO1-L-beta) (EC 1.8.4.-) (Endoplasmic oxidoreductin-1-like protein B) (Oxidoreductin-1-L-beta)
Other Name:	
Human Swiss Prot No.:	Q86YB8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ERO1B
Storage and Stability:	-20°C/1 year